


SPIE Corporate
Member

BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



***Flourescent Label Free 3D CARS Microscopy
for Biomedical Imaging***

Partners:




H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



Company Description

Featured Product: Dual wavelength fs laser system for 3D CARS imaging including tunable Ti:sapphire and Yb fiber laser


Manufacturer of single or double wavelength ultrafast laser systems including ultrashort (ps or fs) pulse, ultrabroadband or broadly tunable Ti:sapphire lasers, Yb-doped fiber lasers, amplifiers and optical parametric oscillators. Their typical applications include time resolved or CARS spectroscopy or nonlinear (2P, SHG or SRS/CARS) microscopy. Manufacturer of ultrafast laser optical coatings including different dispersive mirrors such as chirped mirrors. Complete laser laboratory construction.

H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8



R&D ULTRAFAST LASERS LTD.



FemtoCARS

the

Label-free, 3D Microscopic Imaging System for Real-time in vivo Diagnostics




H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



FemtoCARS

the

Label-free, 3D Microscopic Imaging System for Real-time in vivo Diagnostics






Fig. 3 CARS-images of murine white adipose tissue.



BOOTH NUMBER: 8


R&D ULTRAFAST LASERS LTD.



FemtoCARS


The Label-free,
3D Microscopic Imaging System
for Brain Research

Figure 1. CARS image of myelin fibers in the white matter of rat. In vivo preparation.



BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



FemtoCARS

The concept

BSu3A.28.pdf

Biomedical Optics and 3D Imaging OSA 2012

A simple, cost efficient fiber amplifier wavelength extension unit for broadly tunable, femtosecond pulse Ti-sapphire lasers for CARS microscopy

A. Kolonics^{1,2}, D. Csáti¹, P. Antal¹, R. Szipöcs^{1,2,*}

¹ Research Institute for Solid State Physics and Optics, P.O. Box 49, H-1525 Budapest, Hungary

² R&D Ultrafast Lasers Ltd, P.O. Box 622, H-1529 Budapest, Hungary

* r.szipocs@szipocs.com

Abstract: An inherently synchronized Yb-fiber amplifier based extension unit for femtosecond pulse, broadly tunable Ti-sapphire lasers is introduced, which is well suited for coherent anti-Stokes Raman scattering microscopy.

OCIS codes: (140.7090) Ultrafast lasers; (180.4315) Nonlinear microscopy





Fig. 1 Setup of the CARS extension unit

H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



FemtoCARS

The concept






Fig. 1 The energy diagram of the CARS process: ω_p – the frequency of Ti-sapphire laser (pump), ω_s – the frequency of Yb-amplifier (Stokes), ω_{as} – the frequency of light generated during the CARS process, Ω_{vib} – vibration frequency of the investigated molecule

BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



FemtoCARS

How to make a CARS upgrade for an existing two-photon microscope?

- Add a CARS wavelength extension unit to your femtosecond, tunable Ti:Sapphire laser
- Replace dichroic beam splitter according to excitation and signal wavelengths
- Select laser blocking filter according to excitation wavelength in order to protect detector unit

What will you need?

FemtoCARS Unit



FemtoFiber Yb-fiber amplifier




H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com


BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



FemtoCARS

How to make a CARS upgrade for an existing two-photon microscope?




H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8



R&D ULTRAFAST LASERS LTD.



FemtoCARS

the

Label-free, 3D Microscopic Imaging System for Real-time in vivo Diagnostics



H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com



NIH Public Access

Author Manuscript

Lab Invest. Author manuscript; available in PMC 2013 April 10.

Published in final edited form as:

Lab Invest. 2012 October ; 92(10): 1492–1502. doi:10.1038/labinvest.2012.109.

Multicolored Stain-free Histopathology with Coherent Raman Imaging

Christian W. Freudiger¹, Rolf Pfannl², Daniel A. Orringer^{3,4,5}, Brian G. Saar^{1,†}, Minbiao Ji¹, Qing Zeng^{6,7}, Linda Ottoboni⁸, Wei Ying⁹, Christian Waeber⁹, John. R. Sims⁹, Philip L. De Jager^{8,10,11}, Oren Sagher⁵, Martin A. Philbert¹², Xiaoyin Xu^{6,7}, Santosh Kesari^{13,14}, X. Sunney Xie¹, and Geoffrey S. Young^{6,7,*}

Abstract

Conventional histopathology with hematoxylin & eosin (H&E) has been the gold standard for histopathological diagnosis of a wide range of diseases. However, it is not performed *in vivo* and requires thin tissue sections obtained after tissue biopsy, which carries risk, particularly in the central nervous system. Here we describe the development of an alternative, multicolored way to visualize tissue in real time through the use of coherent Raman imaging (CRI), without the use of dyes. CRI relies on intrinsic chemical contrast based on vibrational properties of molecules and intrinsic optical sectioning by nonlinear excitation. We demonstrate that multi-color images originating from CH₂ and CH₃ vibrations of lipids and protein, as well as two-photon absorption of hemoglobin, can be obtained with subcellular resolution from fresh tissue. These stain-free histopathological images show resolutions similar to those obtained by conventional techniques, but do not require tissue fixation, sectioning or staining of the tissue analyzed.

KLINIKAI ALKALMAZÁSOK





Figure 1. Stain-free histologic imaging with multi-color CRI. (A) Vibrational spectra of the major constituents of tissue: lipids, protein and water. Arrows indicate Raman shifts at which imaging is performed. (B-D) SRS images of a live C2C12 mammalian cell acquired at the CH_2 -stretching vibration at 2845cm^{-1} (B) and CH_3 -stretching vibration at 2940cm^{-1} (C). Multicolor image (D) generated from images (B) and (C) with the green channel (CH_2 image) showing the cell-body and the blue channel (thresholded CH_3 - CH_2 difference image) highlighting the nuclear morphology including a bright nucleolus. (E-H) SRS images of fresh *ex vivo* brain tissue acquired at CH_2 -stretching vibration at 2845cm^{-1} (E), CH_3 -stretching vibration at 2940cm^{-1} (F), and vibrationally off-resonant showing TPA of hemoglobin at a sum frequency of $23,700\text{cm}^{-1}$ (G). Multicolor image (H) generated from images (E-G) with the green channel (CH_2 image) highlighting cytoplasm and myelin sheaths, blue channel (thresholded CH_3 - CH_2 difference image) showing the nuclear morphology, and the red channel (hemoglobin image) highlighting red blood cells. (I) H&E-stained micrograph from the same region in the brain. (J) Same multicolor image as (H) with a different pseudo-color scheme, chosen to mimic the appearance of an H&E-stained micrograph, illustrates the similar image content and appearance of stain-free images and H&E stained sections. Scale bar, $25\ \mu\text{m}$.




Figure 2.

Multicolor stain-free images of various mouse organs (green: CH₂ image; blue: CH₃-CH₂ difference image; red: hemoglobin image) of (A) brain, (B) heart, (C) kidney, (D) liver, (E) lung, (F) muscle, (G) ovary, (H) skin with *stratum corneum* (top) and *stratum basale* (bottom) and (I) spleen. Scale bar, 25 μm.




Figure 3.

Multicolor stain-free images of various brain regions in a wild-type mouse in comparison with paraffin-embedded, H&E and Luxol-stained sections. (green: CH_2 image; blue: CH_3-CH_2 difference image; red: hemoglobin image) of (A) cortex, (B) thalamus, (C) dentate gyrus, (D) cerebellum, and (E) corpus callosum. (F-J) show H&E/luxol stained section of corresponding regions.




Figure 5.

SRS label-free images of demyelination and stroke mouse models. (A) Two-color SRS image (green: CH₂ image; blue: CH₃-CH₂ difference image) of a demyelinating lesion in a mouse with experimental allergic encephalomyelitis (EAE). (A) Perivascular inflammatory cells and histiocytes with large lipid droplets from phagocytized myelin are visible. Red blood cells in a capillary are highlighted by arrows (Hemoglobin image, red, not used in this picture). (B) Two-color SRS image of brain tissue in a mouse stroke model three days post-stroke. Macrophages with many small lipid droplets are dispersed throughout the tissue. Dark areas represent edema fluid. (C,D) H&E stained micrographs of (C) EAE and (D) stroke model from the same region. Scale bar, 25 μm.

Diagnosis of BCC by multiphoton laser tomography

Stefania Seidenari¹, Federica Arginelli¹, Sara Bassoli¹, Jennifer Cautela¹, Anna Maria Cesinaro², Mario Guanti¹, Davide Guardoli¹, Cristina Magnoni¹, Marco Manfredini¹, Giovanni Ponti¹ and Karsten König^{3,4}

¹Department of Dermatology, University of Modena and Reggio Emilia, Modena, Italy,

²Department of Pathology, University of Modena and Reggio Emilia, Modena, Italy,

³Department of Biophotonics and Lasertechnology, Saarland University, Saarbrücken, Germany and ⁴JenLab GmbH, Schillerstrasse 1, 0745, Jena, Germany




Fig. 4. Basal cell carcinoma. (a) 100 μm depth, excitation wavelength 800 nm. Shifting the wavelength to 800 nm, basaloid cells become less visible; employing an excitation wavelength of 820 nm basaloid cells disappear and it is possible to observe empty spaces surrounded by collagen fibres (phantom island); (b) 100 μm depth; (c) 120 μm depth.





Fig. 3. Descriptors of basal cell carcinoma. (a) 50 μm depth with an excitation wavelength of 760 nm. Basaloid cells are disposed in a regular order; they are aligned and elongated in the same direction, showing elongated shape and nuclei, and tightly packed together; (b) 50 μm depth with an excitation wavelength of 760 nm. Sheets of basaloid cells are elongated and aligned in two different directions; (c) 80 μm depth employing an excitation wavelength of 760 nm. At the periphery of basaloid nests, the cells are oriented perpendicularly to the surrounding extracellular matrix (palisading); (d) 90 μm depth with an excitation wavelength of 760 nm. A basaloid nodule (asterisk) is observable as a cell aggregate surrounded by fibres (arrow); (e) 30 μm depth employing an excitation wavelength of 760 nm: sheets of cells intermingled with fibres; (f) 75 μm depth with an excitation wavelength of 800 nm. Increasing the wavelength to 800 nm basaloid cells disappear and only the fibres of the extracellular matrix are visible (phantom island).

Basalioma vizsgálata – 3D FiberScope fejlesztése



Kollagén szerkezet – SHG
Sejtek - autofluoreszcencia

- Lézer hullámhossza
- Színszűrők kiválasztása

Szállézer specifikálása
Optikai szál specifikálása


Biztonságtechnikai vizsgálatok



FiberScope, a kézben tartott nemlineáris mikroszkóp


BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



FemtoRose 100 TUN/NoTouch

The First
Broadly Tunable, femtosecond pulse Ti:sapphire laser



Pioneering Ultrafast LaserTechnology by R&D

INTRODUCING THE FIRST ULTRABROADBAND CHIRPED MIRRORS FOR BROADLY TUNABLE FEMTOSECOND LASERS

526 OPTICS LETTERS / Vol. 22, No. 8 / April 15, 1997

Ultrabroadband chirped mirrors for femtosecond lasers

E. J. Mayer, J. Möbius, A. Euteneuer, and W. W. Rühl

Department of Physics, Philipps University, Riedelstr. 1, D-35032 Marburg, Germany

R. Szipöcs

R&D Laser-Optics Bt., P.O. Box 622, H-1039 Budapest, Hungary

Received November 25, 1996


We report on the performance of widely tunable femtosecond and continuous-wave Ti:sapphire lasers that use a newly developed ultrabroadband mirror set. The mirrors exhibit high reflectivity ($R > 99\%$) and smooth variation of group delay versus frequency over a wavelength range from 600 to 1000 nm. Mode-locked operation with pulse durations of 85 fs was achieved from 600 to 918 nm with only one set of ultrabroadband mirrors. © 1997 Optical Society of America

H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8


R&D ULTRAFAST LASERS LTD.



FemtoRose 100 TUN/NoTouch

the

Broadly Tunable, femtosecond pulse Ti:sapphire laser




- Stable, easy mode-locking (with starter electronics)
- Soliton-like, nearly transform-limited pulses
- Patented Ultrabroadband Chirped Mirror™ optics
 - single optics set from 680 to 1040 nm
- Built in Millennia™ / Verdi™ / Finesse™ pumping (6W, 8W, 10W) – diode-pumped stability
- Sealed, purgeable enclosure
 - reliability, full wavelength coverage
- 15 years of experience
- Labview interface program
- Turn-key, truly hands-off operation (automatic cavity control)

H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



FemtoRose 300 TUN/NoTouch

The Cost Efficient

Long-Cavity, Broadly Tunable, Femtosecond Pulse Ti:sapphire Laser



Pioneering Ultrafast Laser Technology by R&D

INTRODUCING THE FIRST TUNABLE, LONG CAVITY
FEMTOSECOND PULSE Ti:SAPPHIRE LASER

Appl Phys B (2012) 107:17–22
DOI 10.1007/s00340-012-4838-7

Applied Physics B
Lasers and Optics

Tunable, low-repetition-rate, cost-efficient femtosecond
Ti:sapphire laser for nonlinear microscopy


P.O. Antal · R. Szipocs

H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8


R&D ULTRAFAST LASERS LTD.



FemtoRose 300 TUN/NoTouch

The Concept

Schematic of the oscillator




L: pump focusing lens, IC: input coupler mirror, Ti:Sa: titanium-sapphire crystal, BRF: birefringent filter for tuning, P: prisms, HCM: Herriott-cell mirrors, OC: output coupler, S: slit for hard-aperture KLM

Reference

Antal P, Szigligeti A, Kolonics A, Szipöcs R; Tunable, Low Repetition Rate, Femtosecond Pulse Ti:Sapphire Laser for In Vivo Imaging by Nonlinear Microscopy; In: Optics in the Life Sciences Congress (OSA, 4-6 April 2011, Monterey, CA) Paper JTUA12, 2011


Typical measured output power vs. wavelength (at 2.6 W pump)



Two different output couplers were used for short wavelengths (SW OC, crosses) and for longer wavelengths (MW OC, circles).

BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



FemtoRose 300 TUN/NoTouch

The Cost Efficient

Long-Cavity, Broadly Tunable, femtosecond pulse Ti:sapphire laser



Key features

- Low pump laser cost (~ 2.6 W pumping)
- Low, 22 MHz repetition rate
- Higher fluorescence signal
- Lower thermal damage in sample
- No extra-cavity chirp control is required
- Wavelength control by a Zeiss 2P microscope

Applications


- Multiphoton microscopy
- Ultrafast spectroscopy

H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.




FemtoRose 300 TUN/NoTouch

Low Average Power, High Quality Imaging in Two-Photon Microscopy


Tunable, low-repetition-rate, cost-efficient femtosecond Ti:sapphire laser for nonlinear microscopy

Fig. 6 Two-photon absorption fluorescence raw images of mouse dorsal skin using (a) the 76 MHz laser and (b) the 22 MHz laser, at nearly the same excitation power (3.081 mW for the 76 MHz laser and 3.015 mW for the 22 MHz laser). (c) and (d) show the corresponding intensity profiles along the white horizontal line in the middle of the images



BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



FemtoFiber

The

All-Fiber, Femtosecond Pulse Yb-fiber Oscillator/Amplifier System



Pioneering Ultrafast Laser Technology by R&D

INTRODUCING THE FIRST ALL-FIBER, ALL NORMAL DISPERSION, FEMTOSECOND Yb-FIBER LASER

All-fiber, all-normal dispersion ytterbium ring oscillator

J. Fülöp,¹ A. Csörgő,² and B. Kóvács^{1,3}

¹ Research Institute for Solid State Physics and Optics, P.O. Box 49, Budapest 1525, Hungary

² Pazmany Electric Institute of Technology Ltd., 24-26, Követők u., Budapest 1195, Hungary

³ R&D Ultrafast Lasers Ltd., P.O. Box 922, Budapest 1376, Hungary

Received: 2 September 2008; Revised: 9 September 2008; Accepted: 12 September 2008


Published online: 30 September 2008

H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8




R&D ULTRAFAST LASERS LTD.



Pioneering Ultrafast Laser Technology by R&D
INVENTING CHIRPED MIRRORS

Femtosecond Dispersive
and Broadband Optics by IBS technology

Patents




H-1012 Budapest, Attila str. 73, HUNGARY
www.fslasers.com

Phone/Fax: +36 (1) 392 2582
r.szipocs@szipocs.com

BOOTH NUMBER: 8

R&D ULTRAFAST LASERS LTD.



Femtosecond Dispersive and Broadband Optics by IBS technology

- Chirped mirrors (CM)
- Low dispersion ripple, highly dispersive negative dispersion mirrors (MCGTI)
- Ultrabroadband chirped mirrors (UBCM)

